

Z Care

Supervised by:
Dr. Esraa Elhariri
Dr. Asmaa Hashem



TABLE OF CONTENTS

01.

Introduction

02.

Problem Definition

03.

Objective and Motivation

04.

Literature review

05.

Proposed System

06.


Discussion & Experimental Results

07.

Conclusions & Future Work

08.

Our team




Alzheimer

A slowly progressing, but
devastating disease.


01

introduction





INTRODUCTION


Dementia is most common in people over the age of 65. The risk of Alzheimer's disease and other types of dementia increases with age, affecting an estimated 1 in 14 people over the age of 65 and 1 in every 6 people over the age of 80.



INTRODUCTION

Dementia is the name for a group of symptoms associated with an ongoing decline of brain functioning. It can affect memory, thinking skills, and other mental abilities.





55 million


10 million

+

○

○

02. Problem Definition




Problem Definition



Age

lifestyle factors
conditions associated
cardiovascular disease
and with



a family history of the condition

People with Down's syndrome are at a higher risk of developing Alzheimer's disease.

03. Project objective

Project objective



Developing an automated helpful system.

Diagnose Alzheimer's disease in its early stages Accurately using:

Deep learning would give us the ability to

- ※ Slow the disease progression with medicine and exercises.
- ※ Maintain mental function

Make it easier for caregivers to take care of the patient





04.

Literature review

Accuracy 66.7%

Using **MRI**, **FDG-PET**, **CSF**, &
genetic features, with Hybrid (CNN
& RNN).

Weiming Lin.2021



Shangran Qiu.2022

Using **MRI** with CNN.

Accuracy 69.2%



Zhou, P. et al. 2020
PET and MRI images and
distinguish between severity using
CNN and SVM

Accuracy of 87%

Accuracy of 91 %

MRI, PET & CSF data with Using
SAE,DBM & SVM

fnagi.2019



Shen et al.2019

distinguish AD From mild cognitive impairment using PET data Using Deep Belief Network

Accuracy of 86%

Accuracy of 91 %

using **MRI data** with Using autoencoder Accuracy of 91 %

Ju et al.2017



05.

Proposed System



KINEMASTER



Sign Up Functions

9:00

Sign up as a Normal user

Name
Enter your name

User Name
Enter your user name

Phone Number
Enter your password

Email
Enter your email address

Password
Enter your password

Sign up

9:00

Sign up as a Caregiver

Name
Enter your name

User Name
Enter your user name

Patient ID
Enter your user name

Phone Number
Enter your password

Email
Enter your email address

Password
Enter your password

Sign up

9:00

Sign up as a Doctor

Name
Enter your name


User Name
Enter your user name

Phone Number
Enter your password

Email
Enter your email address

Password
Enter your password

Sign up





Login Functions

9:00

Login as a Caregiver

User Name
Enter your user name

Password
Enter your password

Forgot password?

Login

Don't have an account? [Sign up](#)

9:00

Login as a Normal user

User Name
Enter your user name

Password
Enter your password

Forgot password?

Login

Don't have an account? [Sign up](#)

9:00

Login as a Doctor

User Name
Enter your user name

Password
Enter your password

Forgot password?

Login

Don't have an account? [Sign up](#)

9:00

Login Now



Please sign up to continue using our App.

Doctor



Caregiver

Normal User

Password Functions



Languages & Verification code






Normal User Features

The image displays four mobile application screens side-by-side, illustrating various user features:

- Home Screen:** Shows a purple header "Hello, Mohamed". Below it is a "Next Tasks" section with two items: "Drink water" (due at 10:00 AM) and "Lorer amet" (due at 10:00 C). A large circular progress bar indicates "10.00 % Today's progress". At the bottom, there are navigation icons for Home, Note, Plans, and E-call.
- Note Screen:** Shows a purple header "Your Notes". Below it is a "Life" section containing placeholder text about a routine. At the bottom, there are navigation icons for Home, Note, Plans, and E-call.
- Plans & Routine Screen:** Shows a purple header "Plans & Routine". It has tabs for "Plans" and "Routine", with "Routine" being active. Under "Current", there are two boxes: "Everyday Routine" and "Everyday Routine". Under "Old ones", there are two boxes: "Everyday Routine" and "Everyday Routine". At the bottom, there are navigation icons for Home, Note, Plans, and E-call.
- New Plan Screen:** Shows a purple header "New Plan". It includes a "Plan date" section with "Day", "Date", and "Year" options, and a "Tasks" section with a button "Add new task". Below these are two boxes: "Time Plan" and "Time Plan". At the bottom, there are navigation icons for Home, Note, Plans, and E-call.


Patient Features





Care Giver Features

9:00

Hello, Ahmed




Let's go to check out




Add New Patient History



to make Your days productive



Your Medicines



Save Cancel



Today's tasks

- Drink water

Lore ipsum dolor sit amet consectetur...

At 10:00 Am
- Breakfast

Lore ipsum dolor sit amet consectetur...

At 10:30 Am
- Medicine

Lore ipsum dolor sit amet consectetur...

At 10:45 Am
- Reading




Lore ipsum dolor sit amet consectetur...

At 11:15 Am
- Message


Lore ipsum dolor sit amet consectetur...

At 11:15 Am

Emergency Functions




Accounts Profiles



Normal User



Care Giver



Patient

Doctor Features

9:00



Hello, Doctor



search



Search for Patients

3:39



Hello, zcare family



4

Patient Id : 4

Sample ID	Class	Date
-----------	-------	------

No Result

Pass this sample to the machine learning model?

Yes

Cancel



3:39



This operation would take up to 10 mins in Total

3:39

100 KBS WiFi 100

3:44

0.53 KBS WiFi 100

3:45

3.00 KBS WiFi 100



Read Files



Preprocessing

3

Prediction



Read Files



Preprocessing



Prediction



Read Files

Done



Prediction

3:39

3.00
KBS/

Hello, zcare family





4

Patient Id : 4


Sample ID	Class	Date
9	No Result Yet	Sun, 7/9/2023

+ -





-Context Diagram

+ -




-Doctor DFD







-Patient DFD




-CareGiver DFD




-Patient, NormalUser, CareGiver Use case







-Doctor UseCase




Normal User Activity Diagram



-Patient Activity diagram



-Doctor Activity diagram




-CareGiver Activity diagram

+





o

-Sequence diagram



Class diagram




+



o




-ERD Diagram



06.

Discussion & Experimental Results

MRIs



1

data aquisition and preprocessing

2

feature engineering

3

DL classifier

SNPs




feature extraction


feature selection

feature fusion

Framework

Data acquisition preprocessing





Data acquisition

Data preprocessing

GENOTYPE VS PHENOTYPE



GENOTYPE




PHENOTYPE

Difference Between Gene And Allele



2 COPIES





2 ALLELES

↳ CAN BE THE SAME ALLELE



↳ 2 DIFF ALLELES








PURPLE



DOMINANT



GREEN



RECESSIVE

HETEROZYGOUS

↳ FUR COLOUR

↳ ONE ALLELE WILL BE DOMINANT)

↳ OTHER WILL BE RECESSIVE

GENE

ALWAYS

EXPRESSED




Gene	Chromosome	Start	End
PSEN2	1	225,114,896	225,160,427
PSEN1	14	72,662,932	72,766,862
APOE	19	50,090,879	50,114,490
APP	21	26,164,732	26,475,003

filter based on chromosome number

Filter based on position on each chromosome

Feature selection



SNPs



Chromosomes responsible for AD

+



o

APP

PS-1

PS-2

APOE4


4 AD genes in the 4 mentioned chromosomes

Certain genes can increase the risk of developing dementia, including Alzheimer's disease

1. The **APP gene** provides instructions for making a protein called amyloid
2. The **PSEN1 gene** provides instructions for making a protein called presenilin 1
3. The **PSEN2 gene** provides instructions for making a protein called presenilin 2.
Presenilin 2


PSEN2,1 helps process proteins that transmit chemical signals from the cell membrane into the nucleus. Once in the nucleus, these signals turn on (activate) genes that are important for cell growth and maturation

4. **APOE ε4** increases risk for Alzheimer's and is associated with an earlier age of disease onset in certain populations. About 15% to 25% of people have this allele, and 2% to 5% carry two copies **APOE4 is the strongest risk factor gene for Alzheimer's disease** apolipoprotein E gene called APOE4



Feature engineering

Random Forest Classifier



Feature selection using random forest for SNPs


```
[‘rs1150895’, ‘rs12126925’, ‘rs12405469’, ‘rs1295640’,  
‘rs16846644’, ‘rs1783016’,  
‘rs1783025’, ‘rs2014146’, ‘rs2073489’, ‘rs2156079’, ‘rs216762’,  
‘rs2234983’, ‘rs2256331’,  
‘rs2313007’, ‘rs2802268’, ‘rs2829973’, ‘rs2829997’,  
‘rs2830000’, ‘rs2830008’, ‘rs2830012’,  
‘rs2830028’, ‘rs2830033’, ‘rs2830036’, ‘rs2830038’,  
‘rs2830044’, ‘rs2830052’, ‘rs2830088’,  
‘rs362350’, ‘rs362393’, ‘rs373521’, ‘rs3787620’, ‘rs380417’,  
‘rs439401’, ‘rs440666’,  
‘rs452987’, ‘rs466448’, ‘rs466609’, ‘rs7276036’, ‘rs7278838’,  
‘rs8006497’, ‘rs8106922’,  
‘Sex’, ‘Age’]
```

Produce 41 feature


MRIs

preprocessing







CN





AD





MCI







CN



MCI



AD



Model	Accuracy		AUC		Loss	
	Train	Test	Train	Test	Train	Test
CNN 1	0.9661	0.9298	0.9938	0.9810	0.1083	0.2655
CNN Reduced	0.9511	0.9298	0.9912	0.9834	0.1400	0.2014
CNN with Attention	0.9984	0.9561	1.0000	0.9897	0.0101	0.1348

ZCARE Implementation, Used Technologies and Techniques



TensorFlow


colab

+

○

○

07. Conclusions & Future Work



Conclusion

1. Use multimodal data SNPs and MRIs
2. binary classification CNN achieved an accuracy of 89.47% and 91.55% for NC vs. AD and MCI vs. AD, respectively. Finally, for multimodality, ZCARE achieved an accuracy of 79.17 %. Which is high compared to other studies
3. Gives more time for the family and loved ones to understand the disease and plan for their future.



Future work

- 1. Increase SNPs samples**
- 2. Insert personalized medicine with our study.**
- 3. In addition to MRIs and SNPs we intend to use clinical data to improve accuracy**
- 4. Suggesting suitable medicine for each person depending on personalized medicine.**
- 5. Follow the stages of disease development.**
- 6. Add an expert user to manage data and model.**
- 7. Apply Time Series Analysis to understanding the progression of AD for the Patient**

08. OUR Team

Z CARE



- ※ **Yasmeen Saad**
- ※ **Mohamed Ashraf**
- ※ **Shrouk Nasser**
- ※ **Sameul Adel**
- ※ **Mohamed Yasser**

THANKS!

